

**TANNER**  
scientific®

# TN50 & TN50UV CRYOSTATS



WATCH VIDEO!



**IN-STOCK & READY-TO-SHIP!**

## Product Description:

The German-made **Tanner Scientific® TN50 Cryostat** features a smooth microtome for easy, intuitive operation. For extra convenience the system includes a spacious stainless steel chamber, 24 easily accessible freezing positions and a well illuminated cooling chamber that can be effectively cooled down to  $-35^{\circ}\text{C}$ . For additional user safety the system can optionally be equipped with an efficient UV-C disinfection system.



MADE IN GERMANY 

## item # TNR TN50

**Dimensions (WxDxH):** 26 $\frac{3}{4}$ " x 23 $\frac{5}{8}$ " x 43 $\frac{3}{4}$ "

**Weight (w/o accessories):** 280 lbs

**Electrical Parameters:** 115V; 50/60 Hz

- **Exceptional German engineering**
- Ideal for Mohs surgery and traditional pathology applications
- 2 Year Warranty\*
- Optional UV-C disinfection solution
- Removable Flat Mounting Surface Accessory
- Smooth rotary microtome
- Spacious cooling chamber for efficient and safe workflow
- Glass anti-roll plate
- 24 easily accessible freezing positions
- Heated viewing window
- Easy-to-reach electronic control panel
- Automatic defrost
- Removable condensate bottle with level sensor

All trademarks and brands are the property of their respective owners. While every effort has been made to ensure that the material contained herein is accurate at time of publication, Mercedes Scientific is not responsible for typographical or photographic errors, inaccuracies or misprints. Availability and pricing of all products is subject to change or discontinuation at any time without prior notice. Mercedes Scientific reserves the right to change prices, or not sell an item, due to error in pricing.

**TANNER**  
scientific®



TannerScientific.com



888.708.5233

All products are subject to change, availability or discontinuation at any time without prior notice. All trademarks and brands are the property of their respective owners.